



31P-MRS Loader: First time setup

On running “31P_Analysis_Tool.exe” for the first time the following GUI should appear:




To take advantage of the facility to save details of scans and run them in batch mode you must set two environment variables. These are simply “shortcuts” to a specific folder located on a drive accessible to your computer; it is recommended that the drive be permanent i.e. not a removable HDD or USB drive.

To do this, select “Set environment variables” in the “Settings” menu.



You will be prompted to select or create two folders, one which will be assigned to %31PScanCollections%, the other to %31PCollectionResults%.



Upon completing this step, you will be prompted to restart your computer. You may choose to do this at a later point but this must be done before trying to use “Save scan details”, “Save and close” or “Run multiple”.

